
REACTOR UNIVERSAL PARA PLANTA PILOTO
UNIVERSAL REACTOR FOR PILOT PLANT

RE-605
RE-620

1001710
1001711

Indice

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1. General Information

- Handle the parcel with care. Unpack and check that the contents coincide with the packing-list. If any part is damaged or missing, please advise the distributor immediately.
- Do not install or use the equipment without reading this handbook before.
- If you have any doubts, please inquiry by calling to your supplier or J.P. Selecta's technical service.
- Do not use the equipment with liquids which can give off vapours capable of making explosive mixtures.

2. Features

RE glass reactor series of variable speed and frequency which is mainly used for synthetic reaction materials: inject high/low temperature cycle liquid into the intermediate layer of the reaction vessel to effect interlayer materials reactions at the desired temperature; The container can be evacuated according to the needs of the reaction to effect materials reactions in the negative pressure state; The constant pressure funnel can be applied to add materials dropwise and uniformly, while the capacitor could be used to recover some materials in the test by distillation, at a controllable recovery speed.

1. The frames are formed by four stainless steel corners, resistant to corrosion, with structure of change and fixation blockade.
2. Made of high borosilicate glass, the glass pieces are characterized by their excellent resistance to corrosion and to high/low temperatures, as well as their low coefficient of linear expansion to withstand greater thermal shock.
3. The variable speed and frequency stirring provides a high torque and ensures a constant speed of rotation.
4. It is easy to discharge the materials in the reaction vessel, as it is of divisible type that can be detached, without a valve for discharging the dead cavity.
5. With the standard connection mode, disassembling is convenient and efficient.
6. Depending on the need, the container lid can be opened at 3, 4, 5, 6, 7.
7. The vacuum gauge is provided together with the device in order to understand the exact degree of vacuum inside the container.
8. In the vacuum state, the negative pressure could reach 0.095MPa.
9. Air type mixing system, mixing technology bearings.
10. Single, double or triple layer stirring palette, according to the clients selection, to achieve a good synthesis reaction.
11. All the parts in contact with the liquid, made of high borosilicate glass or PTFE material, have good anticorrosive properties.
12. The container body and the cover are held together by a PTFE seal, which can be disassembled conveniently.

3. Technical Data

Note: for 110V 50/60Hz equipments
Heidolph Hei-Torque Core stirrer
Capacity: 25L
Motor power: 105W
Dimensions: 70 x 250 x 180 cm
Rotation speed: 20-2000 rpm

	RE-605	RE-620
Capacity	5 litres	20 litres
Volumen double layer	1,5L	6L
Motor power	90 W	250 W
Maximum torque	132N/cm	200N/cm
Dimensions	130 x 70 x 35cm	180 x 62 x 50cm
Stirring speed	0-250 rpm	0-250 rpm
Power supply	220V/50Hz	220V/50Hz

4. Borosilicate Glassware 3.3

- Chemical properties
- Hydrolytic resistance
- Acid resistance
- Alkali resistance
- Temperature resistance
- Optical

Technical features of working use in standard reactors: RE-605 (part no. 1001710) and RE-620 (part no. 1001711).

Vacuum: -0.9 to +0.5 pressure bars

Working temperature: -30°C to +200°C

PTFE gaskets: -100°C a +250°C (maximum working -20 to +200°C)

LUBRICATION:

High vacuum grease

Range Hg: 10-6 mm

Range: -40°C to 200°C invaluable

This grease does not melt and is been composed for lubrication and sealing of glass and wrenches in high vacuum systems with outstanding stability.

GLASSWARE:

Glass is an organic mixture of metal oxides fused together at high temperatures, which upon cooling, solidifies into the clear, rigid, non-crystalline, versatile material known widely across the globe

Borosilicate glass has a very high resistance to attack from water, acids, salt solutions, halogens and organic solvents. Only hydrofluoric acid, hot concentrated phosphoric acid and strong alkaline solutions cause appreciable corrosion of the glass.

Borosilicate glass has excellent thermal properties at both high and low temperatures. The maximum recommended working temperature for laboratory glassware manufactured is 500°C (for short periods of time only). Special care should be taken at temperatures above 150°C to ensure both heating and cooling is achieved in a slow and uniform manner.

5. Parts Description & Location

Figure (1) STANDARD MODEL



- Pos. 1 - Control box*
- Pos. 2 - Rotating motor*
- Pos. 3 - Motor plate
- Pos. 4 - Constant pressure funnel
- Pos. 5 - Stirring bearing
- Pos. 6 - Loading port solid combination & reactor plug
- Pos. 7 - Universal joint
- Pos. 8 - Circulation connection kit
- Pos. 9 - Connecting flange
- Pos. 10 - Discharge valve
- Pos. 11 - Vacuum meter
- Pos. 12 - Capacitor
- Pos. 13 - Clamping coolant (silicone)
- Pos. 14 - Capacitor dispenser
- 15.- Top flange
- 16.- Central flange
- 17.- Structure
- 18.- Reactor body
- 19.- Rotating pallette
- 20.- Tray gasket
- 21.- Adjusting nut
- 22.- Wheels with brake





** * For standard 220V model. For model 110 / 220V, HEI-TORQUE (page 34),
and according to density-speed.*

6. Installation

1. Open the wooden box and read the manual carefully, check that the equipment comes complete according to the packing list and the label name.
2. Remove the reactor body and frame from the box, with the roller facing down and the wheels blocked.
3. Slowly fill the condenser with water to clean, and place at the installation area.
4. The motor mounting bolts are extended on the motor support, the two-bar stirrer, upto the hole. Depending on the position shown in figure (1), explosion-proof motor or simple motor, the explosion-proof electrical control boxes are attached to the rack.
5. Screw the 4-butterfly nut into the upper tray, upto the release position, remove the adhesive tape, remove the container lid and the sealing gasket from the container.
6. Remove the lock nut, the locking sleeve, and the F4 sealing ring from the top of the container.
7. The body is connected to the container, cover all parts with vacuum silicone grease; cover the stirring blades in the container body with the F4 seal of the upper part of the container; Place the stirring blade on the container lid with a stirring bearing, the cover and the concentric body, the 4 screws in the central tray, the butterfly nuts and the blades, tightly fixed.
8. Leave a distance of more than 20mm between the bottom of the agitation blade and the bottom of the container, place the bearings, the fastening assembly, tighten the stirrer bearing nut on the agitation blade and fix it securely.
9. Connect the rotary motor and the agitation blade with a universal joint and fix it. The start-up should be a flexible rotation, without rim.
10. The feed valves, the sensors, the solid feed port plug, the constant pressure funnel and the condenser dispenser are inserted into the container lid; the capacitor inside the capacitor dispenser, all well fixed respectively. The glass is installed according to the position shown in figure (1).
11. The discharge port and the discharge valve for cleaning, covered with vacuum grease, fixed with a flange. Close them and add water for the test.
12. The vacuum gauge is fixed to the frame. The pumping port of the vacuum table is connected to the pumping hole of the condenser; and the pumping port is connected to the vacuum pump.
13. The cooling water inlet is connected below the condenser loop interface, from the top of the capacitor.
14. The circular import is under the reactor body, while the circular export is in the upper part.
15. Check if the rotation speed is constant after activating the equipment, and adjust the concentricity according to the actual conditions.
16. The temperature control device must be firmly connected to the reaction vessel according to Figure (2), without leakage.

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
Panel de control
Hei-TORQUE Value

Panel de control Hei-
TORQUE Precision

**Coolant double effect DIN 29/32 L: 300mm****Addition funnel 250ml, compensated pressure DIN29 / 32****Blind stopper 29/32 for addition funnel****Stopper 29/32 with outlet 10mm at 90°C****Hei-torque CORE for 5L Part number 1001710 110/220V**

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7. Operation

1. Turn the power switch on and the temperature and the speed will be displayed.
2. The speed is adjusted to 0, click on the Rotation key, slowly adjust the speed, the rotation of the agitation blade, and observe if any anomaly occurs.
3. Pour the materials into the container through solid feed or other inputs.
4. Evacuate the container according to the reaction needs. For the triple reaction vessel, the vacuum pump can be applied to evacuate the outer layer, to achieve the thermal insulation functionality.
5. Set the temperature control device according to the temperature shown, to obtain the desired reaction temperature inside the container.
6. Slowly adjust the knob to achieve the desired rotation speed.
7. Let the water or coolant pass through the condenser according to the process requirements.
8. After the reaction, first turn the heating or cooling temperature control device function off, maintain the cycle and then close the circulation function when it is near to the ambient temperature.
 - Note: The temperature control unit cycle is turned off. The circulation valve can not be closed, leave it open. Otherwise, if the body of the container is too hot or too cold, it could cause inflation and cause damages to the equipment.
9. After the speed test is set to 0, click on the key 1 minute after closing the reactor power.
10. Open the discharge valve to discharge the materials. If necessary, the valve parts or the sealing parts of the container should be cleaned and then reassembled safely.
 - The glass parts must be cleaned after returning to normal temperature. The droplets on the inside and outside should be cleaned after washing to avoid rupture of the container, caused by an irregular temperature or by a sudden change in temperature.

8. Maintenance

The warranty of the equipment against manufacturing problems is of 12 months from the date of purchase. If the equipment has been used improperly, compensation for the damage caused will be charged.

Contact the Technical Service before obtaining professional maintenance to return the equipment to the factory:

1. Clean the machine so as not to damage the maintenance personnel;
2. Provide the description of the equipment failure;
3. Pack the machine properly to avoid damage during transport.

9. Accesories and Spare Parts

Digital thermometer «TC-9226-A».
Temperature range: from -200°C up to 1370°C
Part number 5903030

Criohermostat «FRIGITERM-TFT-30»
Accessory for the coolant circuit
For regulable temperatures from -30 up to 100°C
Part number 3000547

Gable lid
Thermostatic bath accessory which limits evaporation
Part number 3001295

Thermostatic bath DIGITERM-TFT-200 -20
Accessory for the reactor
To regulate temperatures from ambient +5°C up to 200°C
Part number 3000550

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Water recirculator "Interheater"
Accessory of the coolant circuit
To regulate temperature from +5° up to 90°C
Part number 6001491


Rotary vacuum pump
Flow volume 3.6 m³/h. Vacuum limit 0.1 mbar.
Part number 5900620
spare oil for pump
Part number 5515030

Vacuum mebrane pump N-820
Special for solvents
Flow volume 20 litres/minute. Vacuum limit 8 mbar. PTFE rotors. PTFE coated membranes, valves in FFPM, motor protection level IP-44
Part number 5910009


Scrubber unit
Glass trap in Pirex to avoid aspiration of liquids by means of the vacuum pump
Part number 4001611

Universal clamp
In dur-al with jaws in cork coated
Part number 7000403

Double nut
In chrome plated
Part number 7000408


Double glass valve
With three vacuum connection points (coolant, load & drip)
Part number 47078


Glass T
For vacuum connection
Part number 47077


Silicone pipe
External diameter 14 mm. Internal diameter 8 mm
Part number 46041

Corrugated tube
Reinforcement for the silicone tube
Part number 46079

Fixing clamps
Bags of 100 units
Part number 28006 (unit)

PTFE bellows diameter 29/32 x 100mm
Part number 1001825

Conical connection cover 29/32 PTFE
Part number 28515

Glass sensor adapter for reactor
Combines with PH/T°C analyzer part number 4120800
& sensor part number 120801
Part number 4120802

Guide for PTFE stirrers 29/32 Ø 10 MM
High chemical resistance
Low vacuum (5mmHg)
Pressure (3-5 psi)
High performance PTFE/PEEK/glass precision bearing seal
Internal anti-whip sleeve
Part number 1001788

Addition funnel 250ml PTFE key Step 4mm
Part number 1001795

Load faucet 29/32 Step 4mm Olive 12/13
Part number 1001793

Stainless steel funnel with handle 200 mm
Part number 1000909

Support accessories trolley
Total height 90 cm
Lower platform dimensions: 70x70 cm
Upper platform dimensions: 70x40 cm
Part number 1001744 (unit)

Stainless steel bucket 15 litres
Part number 1001600

Large opening Hoffman clamp
Part number 1000064

50 litres reservoir
Part number 0106006

OTHER ACCESSORIES:

- 0046030 PTFE tube 10X8 (1 meter)
- 0046149 Silicon tube to unload 15x20. (1 meter)
- 0046150 Silvar tube 14.4x8 mm reinforced with fiberglass (5 meters)
- 0046156 Silvar tube 19x12 mm reinforced with fiberglass (5 meters)
- 0046084 Viton tube 10x14 mm (1 meter)
- 0046085 Viton tube 12x17 mm (1 meter)
- 0046041 Silicon tube 14x8 mm (5 meters)
- 9515031 Fluid 47 V50 4 litres (silicone oil)
- 1001853 PTFE corrugated tube DN12 connection Ø10 (1 meter)
- 1001915 High vacuum grease 50 gr.

10. Examples of Use of Support Trolleys

Frigiterm TFT & Interheater

Digiterm TFT & Interheater

Digiterm TFT, Vacuum pump,
Scrubber & steel cubes

Vacuum pump, Scrubber
and 50 litres reservoir



11. Problems & Troubleshooting



Failure	Cause	Troubleshooting
There is no power in the equipment	No electricity supply	Check the power source
	Fusion of protective tube	Change the protective tube
	Switch power supply failure	Change the switch
	Circuit board failure	Contact SAT
No display or temperature error	Display panel failure	Replace the display panel
	The display panel cables fall	Fix the connection cables
	Short circuit or open sensor circuit	Check if the seals have loosened: if loose, re-weld. If not, you must replace the sensor
Speed is not shown or there is no rotation	Engine failure	Contact SAT
	Circuit board failure	
	Slow-running rotary switch out of service	Replace the slow-motion switch
	Speed control key out of service	Replace the speed control knob
Abnormal noise	Insensitive rotation of the stirring blades	Reassemble
	Engine failure	Contact SAT
	Debris inside the container	Check and clean them
Pressure leakage and low vacuum degree	Glasses installed and mounted with holes	Reassemble
	Abrasion of the sealing ring	Replace the sealing gasket
	Improper installation of the sealing gasket	Reassemble the gasket
	Aging of the vacuum hose	Replace the vacuum hose



12. Examples of Assembly Schemes


RE-605

Part number 1001710

Vacuum recirculation system
Part number 1001931