Indice

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Index

1 Packing List and Accessories	30
1.1 Package contents	
1.2 Accessories	
2 Safety	31
2.1 Risks to which the operator is exposed	
2.2 Staff gualification	
2.3 Equipment use	32
2.4 Modifications	32
3 Operating Principle	33
4 Installation	34
4.1 Location	
4.2 Power supply connection	
5 Operation	35
5.1 Keyboard description	
5.2 Starting up	35
5.3 Basic Operation	36
6 Analyze Sample	37
6.1 Basic Mode	37
6.2 Quantitative Mode	
6.3 Kinetics Mode	44
7 Utilities	48
7.1 Set Deuterium Lamp	48
7.2 Set W Lamp	49
7.3 Set Printer	49
7.4 Set Time & Date	50
7.5 Get Dark current	50
7.6 Reset WL	50
7.7 Lamp Life Management	51
7.8 Load Default	
7.9 Lamp Change	
/.10 Version of Instrument	
8 Maintenance	52
8.1 Lamps replacing	
9 Guarantee	53
10 Technical Specifications	54
11 Declaración de conformidad CE / Conformity Declaration CE	55

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1 Packing List and Accessories

1.1 Package contents

- Spectrophotometer
- Power connection cable
- Instructions manual
- PC communication software
- Equipment protection cover
- Glass cuvettes
- Quartz cuvettes (2)

1.2 Accessories

- Cuvettes, cuvette holders and supports
- Flow cuvette
- Cuvette holder from 5 to 100 mm
- Thermostatted cuvette holders
- Tube holder

Safety material:

•	Tubes extraction clamps	7000532
•	Heat and spills protection pad glove	1001255
•	Anti-acid gloves	5000044

Safety glasses 5028192

2 Safety

Spectrophotometers UV-2005 include the safety measures suitable for a safety operation in a laboratory.

This equipment has been designed to be used by technicians who have been qualified for its use. For the correct and safety use of this equipment it is important that the laboratory staff know the safety general procedures, together with the safety precautions explained in this manual.

Some of the chemical products used in this spectrophotometer are corrosive and/or inflammable and the samples can be radioactive, toxic or potentially infectious. Be careful to follow the safety standardized procedures in laboratories for samples and chemical products manipulation.

The risky situations are indicated along this manual by means of the safety icons and together with some instructions that must be respected.

2.1 Safety icons

They identify risky situations and the safety measures to be taken.

The icons make reference to the paragraph marked in grey.

Danger risk Risk of danger.

Nok of daliget.

Please, respect the instructions indicated to do the described operation.

Corrosive substances contact risk

There's a contact risk of hard bases or acids that easily attack the skin, cloths and other materials.

Please, respect the instructions indicated to do the described operation.

Electric risk

Risk of electric accident when accessing the areas labeled with this signal or when doing the operations indicated with this icon in the present manual.

Please, respect the instructions indicated to do the described operation.

Burn risk due to contact with high temperature areas The temperature in the area indicated by this icon can exceed 60°C. Please, use thermal resistant gloves to do the described operation.

Please, respect the instructions indicated to do the described operation.

Gloves use compulsory

It is compulsory the use of gloves to do the described operation.

A recommended model is indicated in the «1.2 Accessories» section.

Eyes protection use compulsory

It is compulsory the use of safety glasses or complete facial protection to do the described operation.

A recommended model is indicated in the «1.2 Accessories» section.

Important information

• Important information to obtain good results or for an excellent operation of the equipment.

• Important information to lengthen the equipment life or to avoid the degradation of any of its components.

2.1 Risks to which the operator is exposed

The working process with the described equipment exposes the operator to the following risky situations:

- Corrosive and / or inflammable substances manipulation.
- Sample manipulation which can be radioactive, toxic or potentially infectious.
- Glass pieces manipulation.
- Electric risk.

2.2 Staff qualification

This equipment can only be used by staff qualified in knowing dangers exposed in chemical analysis laboratories.

This equipment can only be used by conveniently trained staff who have read the instruction manual or have been qualified in the use of this equipment.

2.3 Equipment use

This equipment is expected to be used in clinical, pharmaceutical and biochemical laboratories.

The use of this equipment in a way not indicated in this manual can commit the protection guaranteed by the equipment.

This equipment is NOT prepared to work in explosive atmospheres.

2.4 Modifications

The modification of the equipment safety systems operation or manipulation by personnel not authorized by the manufacturer can expose the operator to risks not foreseen in this manual.

3 Operating Principle

UV-2005 equipment is composed of five parts:

- An halogen lamp and a deuterium lamp as light sources.
- A monochromator to select the wavelength and to eliminate second order radiations.
- A compartment for sample cuvettes.
- A photodetector to receive light transmittance and to convert it to electrical sign.
- A display where absorbance, tramittance and concentration values are shown.

Source

chromator

Compartment

Figure 1

The halogen lamp light (or the deuterium one if we work in a wavelength lower than 340 nm), focus on the monochromator entrance where a collimated mirror points the light beam towards the diffraction grating.

The diffraction grating diffuses the light beam to produce a spectrum. A part of it focuses in the monochromator output slot by means of another collimated mirror.

The wavelength light beam selected passes through a filter which helps eliminating second order radiations of the diffraction grating.

This light beam goes through the sample, and then the photodetector produces an electrical sign used to make the transmittance or absorbance calculations, shown in the digital display.

Figure 2

4 Installation

Once the equipment is unpacked, check it is all right. Preserve the packing for some days.

4.1 Location

Place the equipment over a flat, levelled and stable surface, suitable for the equipment weight and not exposed to direct sunlight.

With the aim of having the best equipment performance, keep it isolated, as much as possible, from any magnetic or electrical field source or any other electrical device that could generate high frequency fields.

Install the unit in an area free of dust, corrosive gas and hard vibrations.

Place the equipment at a distance lower than a 1.5m from a power socket.

To work in a comfortable way, left a free space of 50cm around the equipment and a minimum of 20cm behind it.

Remove obstructions or materials that could obstruct the air flow around and underneath the equipment.

4.2 Power supply connection

Please, select a socket near the equipment and suitable for its power.

Use only the supplied cord or similar ones to connect the equipment.

For your safety, the equipment must be earthed.

Before switching on, make sure that the voltage and the power frequency correspond to the indicated in the equipment.

Power point Fuse Switch

USB and printer outlet

5 Operation

The control panel allows the equipment functions selection and informs about the current operating status.

Once plugged to the mains, press the main switch and the equipment will start (ON).

Let it heat up for 20 minutes before taking any reading.

The user can make any operation by pressing the keys.

Any results and the control information will be shown in the display.

5.1 Keyboard description

[START/STOP]	To start or to stop the readings (kinetics).
	To save the readings (other techniques).
[RETURN]	To return to the earlier display or cancel the operation.
[ENTER] To confir	m the data introduced or selected.
[CLEAR]	To clear character or clear the display.
[SET]	To program ABS, %T or energetic levels readings.
[GOTO]	To establish the wavelength
[ZERO]	To make zero ABS (%T)
[PRINT]	To print results.
[0- 9]	Numeric entries.
[.]	Decimal point entries.
[-]	- symbol entry.
[»] [«]	Cursors for selecting and moving through the items.

Keyboard

5.2 Starting up

1. The instrument will position filter first, then initialize lamp position and AD converter. Then D2 / W lamps warm up, Wavelength reset and get dark current. After that will initialize printer, check battery and set up system. And then, the screen will display the main menu.

The self-testing process is as fig 3.

2. If the data in memory has been lost, the instrument will directly calibrate system without any choice for you.

Note: DO NOT OPEN SAMPLE COMPARTMENT LID DURING SELF-TESTING.

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5.3 Basic Operation

Blank

- Push the blank cuvette into the light path.

- Press the key **[ZERO]** for blanking.

Note:

1. If the reference solution is too thick, the energy will be low.

2. Blanking is automatic after a wavelength change.

DO NOT OPEN SAMPLE COMPARTMENT LID DURING BLANKING.

Set Wavelength

(Example: set wavelength in "Basic mode").

- Press [GOTO] (Fig.4).

Use numeric keypad to input wavelength (e.g. 500nm).

- Press **[ENTER]** to change the wavelength, and then blank. After blanking, the screen displays as the upper one (figure 5).

Note:

1. If you don't want to change the wavelength, you can press [RETURN].

2. The wavelength range is 190-1100, other figure is invalid.

Measurement

(Example: set wavelength in "Basic mode").

- Push the blank cuvette into the light path. Then press [ZERO].

- Press **[START/STOP]** in the upper figure on the left to go into testing interface (the lower one of figure 5).

- Press **[START/STOP]** again in the lower one of Fig. 5. The testing result displays on the screen. Repeat this step and the results will be displayed on the screen one by one.

Print test report

(For example: Print the report in "Basic mode", the lower one of Fig. 5).

- Press **[PRINT]** in the lower display of Fig.5 to print the report (Fig. 6).

Note:

After printing, all the data displayed on the screen will disappear.

Before measurement, make a blank reference solution by filling a clean cuvette (or test tube) half full with distilled or de-ionized water or other specified solvent. Wipe the cuvette with tissue to remove the fingerprints and droplets of liquid. Fit the blank cuvette into the 4-cell linear changer and place the cuvette in the slot nearest you. Push the changer so that the cuvette is in the light path (Push the rod in). Close the lid.

6 Analyze Sample

We provide three test methods, which are basic mode test, quantitative mode and kinetics mode.

6.1 Basic Mode

Push the blank cuvette into the light path. In the main menu (the last one of Fig. 3), move the cursor onto "Basic mode" and then press **[ENTER]** to go into basic mode test. After automatically blanking, it will display as the upper one of Fig. 4 and wait for the operator.

Press [RETURN] to exit.

Select Test Mode

You can select between three test modes (T%, Abs, Conc. / factor) by pressing **[SET]** to make choice.

Use [*] and [*] to move the cursor on **Abs.**, then press **[ENTER]**. If the icon " \checkmark " appears at the end of **Abs.**, you have selected successfully. Press **[RETURN]** to exit to the upper screen of Fig. 4, blank at the same time. (Fig. 7)

The operation to select **T%** mode and **Energy** mode, is the same as **Abs** test mode.

Begin to Test

Push the sample into light path, press **[START]** to go into basic mode test screen (Fig. 8).

Press **[START]** again, the test result will be displayed on the screen. The operation is the same if you change your samples.

Note:

1. Five groups of test data can be displayed per screen, a total of 200 groups of data can be saved.

2. In the lower one of Fig. 8, you can change your test wavelength, blank, change test mode and print by pressing: **[GOTO]**, **[ZERO]**, **[SET]** and **[PRINT]** respectively.

You can also clear the test data by pressing [CLEAR].

6.2 Quantitative Mode

It includes two test methods: Standard curve and Coefficient method. In the former method, you can establish a standard curve using standard samples known concentration.

If you know a curve equation, you can select the latter method.

Figure 7

Figure 8

Set Parameters

Press [SET] in the pre-testing interface. It goes into parameters setting interface (Fig.12). You should set Unit, the Number of standard samples, and Concentration of standard samples before testing.

Figure 11

Choose Unit of Standard Samples

Move the cursor onto "**Unit**" by pressing [»] and [«], followed by [ENTER], 8 units are under your selection (Fig 13).

Move the cursor onto the unit you need, then press **[ENTER]**. You can also press the number ahead of the unit directly by means of the numeric keypad. If you don't want to select any item, press **[RETURN]** to go back to the last screen.

Figure 13

Set Quantity of Standard Samples

Move the cursor onto "Number" by pressing [»] and [«], followed by [ENTER]. (Fig. 14)

Input the number of your standard samples by means of the numeric keypad. After pressing **[ENTER]**, it returns to the setting screen.

Note:

Number range is from 1-9, other number is invalid.

Set concentration of standards

Before this step, you must blank first. Pull the blank cuvette into the light path in the pre-testing interface, and press **[ZERO]**. Move the cursor to **Concentration**, and press **[ENTER]**. (Fig. 15)

the numeric keypad, then press [ENTER]. It goes into the next setting interface of

Pull sample 1 in the light path, input the concentration of sample 1 by means of

Figure 14

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January/11

Press **[RETURN]** to go back to pre-testing interface.

Note:

If you want to delete an equation, move the cursor onto it, and just press [CLEAR].

Test and Print

When you have selected the equation, press **[START/STOP]** to go into testing interface in pre-testing interface. (Fig. 19)

After blanking, pull the unknown sample cuvette in the light path, then press **[EN-TER].** The testing result will be displayed on the screen.

Operation is the same if you want to test other samples.

Press [PRINT] to print testing results. (Fig. 20)

the numeric keypad. (Fig.23)

J.P. SELECTA s.a.

Figure 23

Ctra. NII Km 585.1 Abrera 08630 (Barcelona) España e-mail: selecta@jpselecta.es - website: http://www.jpselecta.es January/11

Figure 26

Set Parameters

Press **[SET]** in the coefficient method pre-testing screen, and 3 parameters should be set before testing. (Fig. 24)

Set Coefficient K

Move the cursor on **Coefficient K** by pressing [»] and [«], followed by [ENTER]. (Fig. 25)

Input the new value of **K** by pressing the numeric keypad followed by **[ENTER]**, and it will return to the setting screen.

Note:

The range of **K** is 0 - 9999, other value is invalid.

Set Coefficient B

Move the cursor on **Coefficient B** by pressing [»] and [«], followed by [ENTER]. (Fig. 26)

Input the new value of **B** by pressing the numeric keypad and press **[ENTER]**, and it will return the setting screen.

Note:

The range of **B** is between -9999 - 9999, other value is invalid.

Set Concentration Unit

Move the cursor onto **Conc. Unit** by pressing [»] and [«], and then press [ENTER].

(Fig. 27)

Move the cursor on the unit you need, then press **[ENTER]**. You can also press the number ahead of the unit directly by means of the numeric keypad.

Press [RETURN] to go back to the pre-testing screen.

Blank

Pull the blank solution cuvette into the light path, and then press [ZERO].

Test and Print

Pull the unknown concentration sample into the light path, then press **[START]**, and it will enter the testing screen. Press **[START]** again, you'll get the concentration of the sample. (Fig. 28)

The operation is the same if you have more samples to test.

Press [PRINT] to print the test results. All the data will be deleted after printing.

You can also delete the test results by pressing [CLEAR].

6.3 Kinetics Mode

Move the cursor onto **Kinetics** by pressing [»] and [«] in the main menu, then press [ENTER] to go into kinetics pre-testing screen. (Fig.29)

Figure 30

Set Wavelength

Press **[GOTO]** in the kinetics pre-testing screen, and input the testing wavelength by pressing the numeric keypad (Fig.30).

Press $\cite{SET}\cite{SE$

Figure 32

Set Time Interval

Move the cursor on **Time Interval** by pressing [»] and [«], followed by [ENTER]. (Fig.32)

Input the new interval time by pressing the numeric keypad followed by **[ENTER]**, and it will return the kinetics setting screen.

Note:

- 1. The range of time interval is 1-200, other value is invalid.
- 2. If you don't want to change the value, press [RETURN] directly.

Set Total Time

Move the cursor onto **Total Time** by pressing [»] and [«], followed by [ENTER]. (Fig.33)

Input the new total test time by pressing the numeric keypad followed by **[ENTER]**, and it will return to the kinetics setting screen.

Note:

The range of total test time is 1-120,000, other value is invalid.
 The maximum value of total time is correlative to the interval time, because system only permits that 1000 dots will be adopted at most.

For example:

if the interval time is 1s, then the maximum total test time is: 1 X 1000=1000s; if the interval time is 20s, the maximum total test time is:

 $20 \times 1000 = 20,000s;$

if the interval time is 200s, the maximum total test time is **NOT**: 200 X 1000=200,000s

However, due to the reason that the range of total time is 1-120,000, when the interval time is 200s, the maximum test time will be 120,000s!

3. If you don't want to change the value, press [RETURN] directly.

Set Test Mode

Move the cursor onto **Test Mode** by pressing [»] and [«] in the setting screen, followed by [ENTER]. (Fig. 34)

Move the cursor onto the test mode by pressing [»] and [«], and then press [ENTER] to select. Press [RETURN] to go back to the kinetics setting screen.

Set Upper Limit

Move the cursor onto **Upper Limit** by pressing [»] and [«] in the setting screen, and then press **[ENTER]**. (Fig. 35)

Input the new upper value by pressing the numeric keypad followed by **[ENTER]**, it will return to the kinetics setting screen.

If you don't want to change the value being set last time, press **[ENTER]** or **[RE-TURN]** directly.

Note:

The range or A is -0.3 - 3, the range or T is 0 - 200, other value is invalid.

Set Lower Limit

Move the cursor onto Lower Limit by pressing [»] and [«] in the setting screen,

January/11

Figure 38

followed by [ENTER] (Fig.36)

Input the new lower value by pressing the numeric keypad and then press **[ENTER]**, and it will return to the kinetics setting screen.

If you don't want to change the value being set last time, press **[ENTER]** or **[RE-TURN]** directly.

Note:

The range or A is -0.3 - 3, the range or T is 0 - 200, other value is invalid.

Sample Test

After you have finished setting **Time Interval**, **Total Test Time**, **Test Mode**, **Upper Limit** and **Lower Limit**, press **[RETURN]** to go back to the kinetics pre-testing screen. (The upper display of Fig. 37)

Pull the blank cuvette in the light path, then press **[ZERO]** for blanking. Pull the sample cuvette in the light path, and then press **[START]** to go into the testing screen.

Press [START] again to begin your test. (Fig. 37)

When the testing begins, the curve will be drawn on the screen in real time. And the time will also be displayed in the middle of the bottom screen in real time.

Press **[STOP]** to interrupt the test. When **[START]** is being pressed again, the test will start again from the very beginning.

Note:

The curve cannot be printed, if you want to check the data of every dot, you can choose **Display Data** in the parameter setting screen.

Display Data

Move the cursor onto **Display Data** by pressing [»] and [«] in the setting screen, followed by [ENTER]. (Fig. 38)

Five lines of data will be displayed on every screen. Press [»] and [«] to scroll.

Data Print

Press [PRINT] in the down display of Fig.38, and the data will be printed.

After printing, the data will be deleted from the screen and RAM.

You can also delete the data if you press [CLEAR]. (Fig. 39)

7 Utilities

From the Utilities menu, we can browse and program different parameters. First, move the cursor to **Utility** by pressing [»] and [«] in the main menu and then press [ENTER]. You could set ten items (Fig. 40)

Figure 41

7.1 Set Deuterium Lamp

Move the cursor onto **D2 Lamp on/off** by pressing [»] and [«] in utility setting screen, and then press [ENTER] to go into D2 controlling interface. (Fig.41)

Move the cursor onto **OFF** to shut off the D2 lamp by pressing [»] and [«], and then press **[ENTER]**. Press **[RETURN]** to exit.

Move the cursor onto **On** to lighten the D2 lamp by pressing **[*]** and **[*]**, and then press **[ENTER]**. Press **[RETURN]** to exit.

Note:

The D2 lamp needs about 20 seconds to warm up before being lightened! (Fig. 42) Shutting off D2 lamp is strongly recommended when you don't use it after self-testing.

Figure 48

In the time and date setting screen, press [»] and [«] to move the cursor. You can set them by pressing the numeric keypad followed by [ENTER].

The range of year is 0 - 99, month is 0 - 12, day is 1 - 31, hour is 0 - 24, minute is 0 - 59, week is 1 - 7, other number is invalid.

7.5 Get Dark current

When the ambient circumstances changed, those such as temperature, voltage or environment lightness, you should get dark current before testing.

Move the cursor onto **Dark Current** by pressing [»] and [«] in utility setting screen, and then press [ENTER] to get dark current. (Fig. 47)

It will return to the utility setting screen when the dark current gets finished.

7.6 Reset WL

Please, perform this operation if you suspect that wavelengths selected by the user do not coincide with the ones selected by the equipment in the monochromator.

Move the cursor onto **Reset WL.** by pressing [»] and [«] in the utility setting screen, and then press [ENTER] to reset wavelength. (Fig. 48)

The calibrating time is about 1.5 minutes.

7.7 Lamp Life Management

Move the cursor onto **Lamp Life** by pressing [»] and [«] in the utility setting screen, and then press **[ENTER]**. (Fig. 49)

Move the cursor onto **Display Lamp Life** by pressing [*] and [*] in the lamp life setting screen, and then press [**ENTER**]. You'll find the used time of D2 and W

January/11

(It can be modified without notice)

lamp. (Fig. 49)

When the used time of D2 lamp is higher than 1000 hours, replacing a new one is strongly recommended.

When the used time of W lamp is higher than 1000 hours, replacing a new lamp is strongly recommended.

When you have replaced a new lamp, you'll have to go into the lamp life setting screen, and move the cursor on the "Reset W Life" or "Reset D2 Life" by pressing [»] and [«] followed by [ENTER] to reset them into zero.

7.8 Load Default

Move the cursor onto **Load Default** by pressing [»] and [«] in utility setting screen, and then press **[ENTER]**. (Fig. 50)

If you want to default the parameters, choose "Yes" in the lower screen of Fig. 50. If not, you should choose "No".

Press **[RETURN]** to exit without any choice.

7.9 Lamp Change

Move the cursor on "**Lamp Change**" by pressing [»] and [«] in utility setting screen, then press [ENTER]. (Fig. 51)

Input the new wavelength you want by pressing the numeric keypad followed by **[ENTER]**.

Lamp change wavelength means the point at which W lamp and D2 lamp exchange. The range of it is 300 - 400, other value is invalid.

7.10 Version of Instrument

Move the cursor onto **Version** by pressing [»] and [«] in the utility setting screen, and then press [ENTER]. Then you will find the detailed messages of the instrument.

8 Maintenance

8.1 Lamps replacing

Tungsten-halogen lamp

Switch off and unplug the equipment. Remove the four screws from the spectrophotometer sides. Remove the **sample selection bar** unscrewing counter clock.

Remove the equipment cover very carefully and place it to the right.

BE AWARE of not PULLING OUT the cable.

Unscrew the two lamp chamber screws and remove the cover.

Attention The cover could be very hot, remember to use gloves!

Unplug and remove the lamp from the ceramic base (black connector). Insert the new lamp, pushing as far as possible.

CAUTION:

1) Do not touch the lamp with bare fingers. Use a cloth to handle it.

2) There's no polarity difference between both tungsten lamp pins.

Connect the equipment and adjust the lamp in order that the light beam focuses on the monochromator inlet slot.

Unplug the equipment, cover again the lamp chamber and install the cover. Pay attention not to pinch the cables during the process.

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Attention Use UV protection glasses while changing the deuterium lamp!

Switch off and unplug the equipment. Remove the four screws from the spectrophotometer sides. Remove the **sample selection bar** unscrewing counter clock.

Remove the cover carefully and place it to the right of the equipment.

BE AWARE of not pulling out the PANEL with the cable.

Unscrew the two lamp chamber screws and remove the cover.

Connector

Fig. I

Fig. II

Attention

The cover could be very hot, remember to use gloves!

Unplug the 3 wire connector by pulling it up.

Loosen and remove the two white screws (Fig. I and II in the next image) from the deuterium lamp.

Replace for a new pre-aligned lamp provided by **SELECTA** or an authorized technical service.

Make sure that the lamp bridle wire coincides with the pintle position. Then, tighten both screws.

Switch again the cable connector (make sure that the cable connection orientation is the correct one.)

Switch the equipment on. After self-diagnosis, select 300 nm and press [ZERO]. Check that the light beam focuses the monochromator inlet slot.

As the lamp is a pre-aligned one, it will need just a minimum adjustment. Unplug the equipment, cove the lamp chamber again and tighten both screws. Install again the equipment cover. Pay attention not to pinch the cables during the process.

8.2 Battery replacement

UV-2005 can save many data in the RAM memory, and so a battery is needed in the circuit board.

You will have to replace the battery when it is exhausted.

Switch off and unplug the equipment (VERY IMPORTANT: high voltage).

Unscrew the board's 13 screws and remove it. Remove the old battery and place the new one (Fig. III) Place again the board and the 13 screws.

9 Guarantee

This product is guaranteed for one year from the date of purchase and it covers defective components and assembly.

Please keep the purchase ticket as it is the only valid document for the purchase date.

This guarantee does not cover damage caused by an incorrect use of the equipment. This equipment must be properly validated before its use in the value sample analysis for the user.

Any return without the authorization of J.P. SELECTA, s.a. is not accepted.

10 Technical Specifications

Mains voltage:	230V - 110V 50/60Hz (See features label)
Electric power:	110W (See features label)
Maximum intensity:	0.4 A
Fuses:	3.15 A fast. dimensions: 20x5 mm
Network connection cable:	Phase, Neutral, Earth (section 1mm ²)
Weight (Net):	14 Kg
Dimensions: Depth x Width x Height	370 x 470 x 180 mm
Model	UV 2005
Wavelength range	190-1100nm
Spectral wavelength	2nm
Optical system	one beam with 1200 lines/mm net
Wavelength accuracy	± 0,5nm
Wavelength reproducibility	0,3nm
Wavelength resolution	±0,1nm
Accuracy	±0,3%
Repeatability	±0,2%
Reading range	from -0,3 to 3A,T 0-200%
Diffuse light	0,05% @ 220nm
Stability	±0,002 A / h @ 500nm
LCD screen	(128x64 dots)
Membrane keyboard	22 keys
Reading modes	transmission - absorption - concentration
Photodetector	Silicon photodiode
Sample compartment	Standard cuvettes 10 mm.
Cuvettes capacity	100 mm.
Light source	Tungsten and deuterium lamps
Ports outlet	USB- parallel port (printer)
Sound level	Less than 60 dBA
Environmental conditions:	Inner use
	Altitude up to 2000m
	Temperature between 5°C & 40°C.
	Maximum relative humidity 80% for tem- peratures up to 31°C, linearly decreasing till 50% relative humidity to 40°C
Overvoltage	Category II
Pollution rank	2

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Electrical material for measurement, control and use in laboratory. Electromagnetical compatibility requirements (CEM)

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Electronical equipment safety requirements for measurement, control and use in laboratory. Part 1: General requirements.

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Ramón Ramón Rble. Técnico

David Pecanins Rble. Calidad