MANUAL DE INSTRUCCIONES

Indice

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Index

1 Packing list and accessories	13
1.1 Packing list	13
1.2 Accessories	13
2 Safety	14
3 Overview.	14
4 Technical specification	14
5 Description	15
5.1 General	15
5.2 Rear panel	15
5.3 Using the touch screen	15
6 Installation	15
6.1 Locate and connection to mains	15
6.2 Electrode and temperature probe connection	15
6.3 Measuring temperature. ATC correction	16
7 Operation	16
7 Operation	16 16
7 Operation. 7.1 Initial screen	16 16 16
7 Operation	16 16 16 17
7 Operation	16 16 17 17
7 Operation	16 16 17 17 17
7 Operation	16 16 17 17 17 18 18
7 Operation	16 16 16 17 17 18 18 19
7 Operation	16 16 17 17 18 18 19 19
7 Operation	16 16 17 17 18 18 19 19 19
7 Operation	16 16 17 17 18 18 19 19 19 19
7 Operation	16 16 17 17 18 18 19 19 19 19 19
7 Operation	16 16 17 17 17 17 19 19 19 19 19 19 19

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English Version

1 Packing list and accessories

First, check that the contend is complete and looks in a good state. Keep the transportation box for some days.

1.1 Packing list

•	pH-meter		(1)
•	Arm with probe holder	(2)	
•	Dust cover	(3)	
•	Short-circuited connector.	(4)	
•	Instruction manual	(5)	
•	RS-232 communication cable	(6)	
•	Power cord	(7)	
•	Printer Communication cable	(9)	

1.2 Accessories pH combined electrodes:

Fig.	Code	Range	Range	Temp.	Electrolite	Size	Material	Application
		pН	(°C)	Probe			Body	
1 4	120102	0 - 14	0-80°	No	Liquid	Ø12 x 120	Glass	Laboratory general use.
2 4'	120105	0 - 14	0-60°	Si	Gel	Ø12 x 120	Ероху	Anti-chock body. Laboratory general use and outdoors.
3 4	120101	0 - 14	0-60°	No	Gel	Ø12 x 120	Ероху	Anti-chock body. Laboratory general use and outdoors.
4 4 ⁻	120114	0 - 10	0-60°	No	Gel	Ø12 x 120	Ероху	Surface sensor. Skin, meat, wet semisolids, culture.
5 4 ⁻	120104	0 - 14	0-60°	No	Gel	Ø12 x 120	Ероху	Sharper sensor. Semisolid food, fruits, meat.
6 4 ⁻	120103	0 - 14	0-80°	No	Gel	Ø8x 60	Glass	For in tube measures. Aqueous and solutions.
7 4'	120106	0 - 14	0-80°	Si	Liquid	Ø12 x 120	Glass	Fast response on aqueous means.

2 Safety

Across this manual, some risk of danger situations are remarked. Proceed as described to avoid any hazard.

• Warning symbol of a dangerous operation: Proceed according indicated instructions.

• Important Information symbol: Information important to obtain good results or for an optimal operation of the equipment or to extend the life of the equipment or to avoid the degradation of some of its components.

• Using a pH-meter, the operator, could be exposed to chemical hazard.

• This equipment must only operated by personnel trained a suitable qualified to work in a chemical laboratory and its common hazards and who have read and understood this manual or have been qualified on its manipulation and operation.

• This equipment is not designed and not manufactured to work in an explosive ambient area.

3 Overview

The pH-2006 is a table-top pH-meter with a 320x240 graphic and operated by a touch screen. Suitable for indoor measures on the laboratory.

The pH-2006 also measures temperature for ATC. If no temperature probe is connected, sample temperature can be entered by using the keypad and also ATC is applied on the measure. (ATC: Automatic temperature compensation)

The reading in can be used to take measures with Ion selective and Redox probes. Combined or reference electrodes can be used.

Power supply:	230V 50Hz
pH range:	-218.00
pH precision:	+ 0.01 pH
pH resolution:	0.001 pH
mV range:	±1999.9 mV
mV precision:	±1 mV
mV resolution:	0.1 mV
Temperature correction range:	-5 105 °C
Temperature precision:	±0.4°C
Temperature resolution:	0.1°C
Input impedance:	5 x 10 ¹² Ohm
Calibration points:	2 to 5 points
Buffer recognition: (25°C)	1.68 3.68 3.78 4.00 6.86
	7.00 7.41 9.18 10.01 12.45
Storing capacity:	Up to 200 meassures
Ambient working conditions:	15 30°C humidity (Hr) lower than 70%
Dimensions:	300 x 200 x 70 mm
Weight:	1.5Kg
Output:	RS-232 serial output

4 Technical specification

10

Do not press so much the pointer.

Remove the electrode protector.

Keep placed the electrode protector until their usage.

5 Description

5.1 General

- 1. Electrode.
- 2. Temperature probe.
- Probe holder.
 Graphic touch screen.
- 4. Graphic touch screen.

5.2 Rear panel

- 5. Reference electrode socket (Not combined electrode)
- 6. Electrode socket pH o mV. (Combined electrode)
- 7. Printer/RS-232 connection socket.
- 8. Main switch.
- 9. Earth socket.
- 10. Temperature probe socket.
- 11. Fuse (0.5A).
- 12. Mains socket.

5.3 Using the touch screen

All functions on the pH-2006 are selected by pointing on the touch screen.

To point any option shown on the screen you can do it by your fingers or, better, by the pointer.

Please, in order to not scratch the screen, press the pointer gently.

It is not necessary to press the pointer to much. The screen is enough sensitive to act at gently pointing.

In order to protect the touch screen against excessive dust, keep the pH-meter with its anti-dust cover.

6 Installation

6.1 Locate and connection to mains.

Place the pH-meter on a flat surface, clean of dust, not to direct sun light exposition, and near to a mains plug. For an ergonomic use, let 30cm free space round the pH-meter.

Use the power cord to connect pH-meter to mains plug provided with an earth conductor.

6.2 Electrode and temperature probe connection

The electrode must be stored with its immersion cap filled with electrolyte. The electrode and the temperature probe must be connected on the rear panel.

To connect the electrode, remove the short-circuit connector from the pH-meter and to connect the electrode to BNC socket. Please keep the BNC socket with the short-circuit connector when the pH electrode is not connected on the socket.

Note: The BNC connectors are extracted and placed turning 1/4 turn.

Keep the electrode with its immersion cap until use.

CAUTION:

The immersion cap contend KCI 3M as electrolyte which is corrosive.

Take suitable cautions to avoid electrolyte be in contact with skin or wear. Avoid any contact with eyes. In case to be in contact clean with abundant water.

6.3 Measuring temperature. ATC correction

The pH-2006 includes ATC (Automatic temperature correction). This means that the reading is corrected if the measuring temperature is different from 25°C.

The measuring temperature is the one the temperature probe reads or if the temperature probe is NOT connected it is the temperature set by the keypad.

See 7.2 paragraph to set the measuring temperature by the touch screen: (Temperature probe NOT connected on the rear panel socket)

Note: If the measuring temperature set has a value far from the actual sample temperature, the pH reading could have considerable error.

Keep placed the short-circuited stopper when the electrode is not connected on the pH-meter.

7 Operation

7.1 Initial screen

The initial screen is the pH measuring screen. As well as the pH more information is done:

- System time.
- Measuring temperature.
- mV
- Temperature probe type.
- ATC status.
- Last Buffers used for calibration.

Buttons: «SETUP», «VIEW», «STD», «SAVE», «PRINT», «DEL»

7.2 SETUP screen

Use this menu to set pH-meter general parameters as:

- Current time:
 - After this line is pointed the "Set current time" window is opened:
 - Point the item to change: year, month,... (SETUP button

appears)

- Point several times over the item to change to set it.
- Point «SETUP» to validate the new setting.
- Point on the X on the Top right corner to close the window.
- Resolution (mV):

On pointing over this line option, mV resolution is changing successively.

• Resolution (pH):

On pointing over this line option, pH resolution is changing successively.

• Default temperature.

This will be the measuring temperature when no temperature sensor will be connected.

After this line is pointed, the "Set default temperature" window is opened:

- Point «SETUP» to change de value.
- Point «ENTER» to validate.
- Temp. electrode

After this line is pointed, a new window is opened. Select the type of temperature probe connected on the "TEMP" socket.

Zero (mV) After this line is pointed a new window is opened to adjust the 0mV reference. Usually run this operation connecting the "Short-circuited" plug on the "MEAS" socket.

7.3 VIEW screen

Use this menu to view:

- Last Calibration details.
- Stored readings.
- «Last STD Datum» After this line is pointed the «Last STD Datum» window is opened: It shows all information about the Buffers used for calibration.
 - «¬ ® » Move to next standard Buffer details.
 - «STD» Jumps to Calibration procedure. See XXX
 - «SETUP» Jumps to standard Buffers selection screen.

Point over $\mathbf{X} \mathbf{\hat{p}}$ exit the window and return to initial screen.

«Stored Datum»
 This screen shows all previous reading stored by the user on the pH-meter.

Function Keys:

• «PRINT» If printer (See 7.6) is attached to the pH-meter starts to print all or just the selected reading.

- «DEL» Deletes all stored reading or just the selected one.
- «PgUp» «PgDn» Scroll page.

Point over 🖾 o exit the window and return to initial screen.

7.4 STD screen (pH-Calibration)

Calibration is the process to calculate the relationship between the mV readings and the actual buffers pH value. It is necessary, because all electrodes have different features.

In a few words, calibration steps are:

• Select the "Buffer Group" that is going to use. This means that you must set what Buffers are going to be used for calibration.

• Measure each Buffer with the electrode. On this step, the pH-meter calculates the correction to apply from the pairs measured mV - pH.

From the starting menu, point over «STD» button to proceed with pH-meter calibration.

Note: Usually a pH-meter must be calibrated:

- First time use
- Daily
- After change the electrode
- After measuring strong acids (pH<2)
- After measuring strong alkalis (pH>12)
- After measure very dense solutions that can obstruct the membrane.

7.4.1 Calibration points

The pH-2006 can be calibrated up to 5 points. However 3 points calibration is wide used.

As a suggestion choose the calibration point following this criteria:

- One point is always pH=6.98
- If you are going to measure an acid solution choose pH=4.00
- If you are going to measure an alkali solution choose pH=9.18
- If you are going to measure in the whole 0...14 range use 5 points.
- Choose calibration points, at less, 1 pH unit different between them.

pH-2006 recognizes, automatically 10 buffers. Buffers can have a value different that the ones showed on the "Set buffer group" screen. Because the exact value of the buffer can be entered manually.

To select the "Buffer group" (The buffers that are going to be used on calibration. Point over «SETUP» on the screen «Standarization». A new window appears with the standard buffers. Point 2 times over the desired buffer icon to joint or remove from the group.

7.4.2 Calibration procedure

Note: To get an accurate calibration, be sure that buffers are in a well status (not contaminated, caducity date, temperature, etc...).

1. Adjust Zero mV.: Place the "short-circuit stopper" (pg.13) in the "electrode socket". See rear panel in 5.2 (pg.4). Press "SETUP", select "Zero(mV)" and press "ENTER" to enter mV calibration panel. The meter displays "Adjust zero(mV)?". Press "ENTER" to confirm or "CANCEL" to exit to starting working state.

2. Check, in the bottom "starting screen" that the Buffer group is really what you are going to use. See 7.4.1 to change it.

- 3. On the "starting screen" press «STD».
- 4. Press «Enter» to validate the process.
- 5. Immerse the electrode cap on the first buffer. Wait until the reading becomes stable.
- 6. Press «ENTER».

• The pH-2006, must recognise the Buffer from the "Buffer group", and adds the buffer icon to the bottom window "STD Result".

• The window "End standardization?" is shown. Point over "Continue" or over "End STD" to finish the calibration.

- 7. Clean and rinse the electrode cap with clean water.
- 8. Repeat the steps: 5. to 7. as necessary to calibrate in all desired points.

9. On the window "End standardization?". Point over "End STD" to finish the calibration.

Non standard buffers or buffers at different temperature than 25°C:

Adjusting the entry for Buffers that its value does not match with the actual value:

On the step: 6. instead to press «ENTER», press «Auto» and the «ENTER». The "Type Buffer pH" is shown to, manually, enter the Buffer value and validate pointing over «ENTER».

Attention

Buffer solutions are irritant chemicals. Avoid any contact with eyes and use a lab coat to protect your wear.

In case of chemicals be in contact with skin, clean, quickly, with abundant water.

To prepare the buffer solutions just take a 300-400 ml flask, and solve the pounder bag into 250ml of distilled water. Stirring until solved.

Keep these solutions in a closed flask. Do not use this buffer for calibration after 1 moth after preparation.

7.5 SAVE screen

Save the current result. The measuring data is stored:

- Date and time.
 pH reading.
- Temperature. mV.
- Slope.

7.6 PRINT screen

Sends the reading or the stored data to the printer or to the computer via RS-232 serial communication.

Use printer 4120113 or other serial printer with the following configuration:

9600 b.p.s. No parity.8 bits data 1 bit start 1 bit stop.

Point over «Print current result» or «Print all» to print only the current measure or all the stored measures.

To download data to the computer ask J.P. SELECTA,s.a. for the communication software. It is available and free.

7.7 DEL screen

Deletes all stored data. Use «VIEW» See 7.3 to delete individual stored data.

Attention

After data or datas are deleted, there is no means to recover it.

8 Maintenance

Keep clean and ${\rm d} r {\rm y}$ the electrode socket. Use the short-circuited plug on the socket instead of let it on the air.

8.1 Electrode handling

- If the electrode comes with its own instruction sheet, follow it.
- Keep the electrode with its end cap filled with KCI 3M when not used.

• Take suitable cautions to avoid electrolyte KCI 3M be in contact with skin or wear. Avoid any contact with eyes.

• Avoid than the electrode left more than 1h without the end cap and becomes dry on the air. It can lost its features if becomes dry. For a short period, left it immersed on a glass of water.

· Clean and rinse with abundant water after measure.

• Avoid that the electrode will be immersed in distilled water, protein solution, acidity fluoride solution or silicon oil for a long time. This can cause the membrane obturation and its response will be slowliness and inaccurate. Clean with HCl 10% for 1h to clean it.

• Be sure that mains plug is grounded.

9 Warranty

The warranty covers the defective components or faulty assembling of the unit for 1 year since the purchasing date.

Keep purchasing bill as the only valid document for the purchasing date.

The warranty does not cover the damages caused for a unit malfunction. The unit must be suitably qualified before to be used for valuable samples processing.

Return good are not accepted without previous J.P. SELECTA, s.a. approval.

10 Statement of conformity CE

The manufacturer:

J.P. SELECTA, s.a. Ctra. NII Km 585.1 08760 ABRERA (BARCELONA) ESPAÑA

Declares that the laboratory equipment:

pH-meter pH-2006 Code: 4120600

Meet the following European Directives:

73/23/CEE Directive for electrical safety. 89/336/CEE Directive for electromagnetic compatibility(EMC).

Meet the following Standards:

EN 61326:2002 Electrical equipment for measurement, control and laboratory use. EMC Requirements.

EN 61010-1:2002 Safety requirements for electrical equipment for measurement, control and laboratory use. Part 1: General requirements.

Abrera, June 2006

Ramón Ramón Dtor. Técnico.

David Pecanins Rble. Calidad