



## Incubation chamber "Boxcult"

FAN CONVECTION.

ADJUSTABLE TEMPERATURES FROM AMBIENT +5 °C UP TO 57 °C.

STABILITY:  $\pm 0.25$  °C, UP TO 37 °C HOMOGENEITY:  $\pm 1$  °C, UP TO 37 °C.

SET ERROR:  $\pm 2$  % OF THE WORKING TEMPERATURE. RESOLUTION 0.1 °C.

### SAFETY:

SAFETY STANDARD EN 61010. OVER TEMPERATURE SAFETY THERMOSTAT FITTED.

### FEATURE

Made of transparent methacrylate that allows the user to see inside the incubator during operation. To facilitate the access to the working area the unit has a wide front door, and a removable base made of AISI 304 stainless steel. The fan convection circulation system ensures an even and rapid recovery of temperature.

A 30 mm Ø port at the rear can be used to connect power to apparatus inside the chamber.

Supplied as accessories, the removable base allows the Boxcult to be mounted on the "Rotabit" reciprocal / orbital shaker. (described in the stirrer section.)

The metallic top of the chamber includes the heating elements, air circulation fan and temperature control.

### CONTROL PANEL

Main switch.

Digital electronic temperature control.



*Incubation chamber "Boxcult" Part No. 3000957 with base Part No. 3001172 and support rack with two shelves Part No. 1000973. Supplied as accessories.*

### MODEL

Part No.	Capacity litres	Height / Width / Depth (interior) cm			Height / Width / Depth (exterior) cm			Power W	Weight Kg
<b>3000957</b>	<b>110</b>	<b>50</b>	<b>47</b>	<b>47</b>	<b>61</b>	<b>51</b>	<b>51</b>	<b>220</b>	<b>18</b>

Supplied without bottom base, or stainless steel rack and shelves.



*Orbital & shaker stirrer "Rotabit" part number 3000974 with incubation chamber "Boxcult"*

### ACCESSORIES

**Removable bottom** base made of AISI 304 stainless steel. Part No. **3001172**

**Stainless steel rack** with 4 shelves positions, each one separated by 9 cm. Comes complete with 2 removable shelves. Useful dim. 43 cm long and 41 cm wide. Part No. **1000973**

For an easy handling, all control devices are outside the chamber enclosure