



# Water and oil baths "Precisdig" with or without anti-evaporation folding lid

DIGITAL CONTROL AND DISPLAY OF TEMPERATURE AND TIME.

ADJUSTABLE TEMPERATURE RANGE FROM AMBIENT+5 °C TO 99.9 °C, STABILITY  $\pm 0.1$  °C, HOMOGENEITY  $\pm 1$  °C.

ADJUSTABLE TEMPERATURE RANGE FROM AMBIENT+5 °C TO 200 °C, STABILITY  $\pm 0.5$  °C, HOMOGENEITY  $\pm 2.5$  °C.

**SAFETY:** DIN 12877 STANDARD. THE UNIT WILL CUT OUT IF THERE IS INSUFFICIENT LIQUID OR IF THE BATH TEMPERATURE EXCEEDS THE SET TEMPERATURE. A MANUAL RESET ADJUSTABLE TEMPERATURE SAFETY THERMOSTAT IS PLACED AT THE BACK OF THE UNIT.

**Programmable temperature and time. Interior and exterior made from Stainless steel**

NEW  
DESIGN



## FEATURES

Double skinned tank of AISI 304 stainless steel with an external case of AISI 304 stainless steel.

Heating elements embedded in INCOLOY heat and corrosion resistant stainless steel.

Comes complete with drain tap.



## CONTROL PANEL

General switch with signal lamp.

1. Heating indicators; open thermostat; failure.
2. Reading display (temperature or time).
3. Setpoint value display (temperature or time).
4. Visualization mode change.
5. Increase value.
6. Decrease value.
7. Start / stop button.

## MAKROLON® FOLDING LID PROPERTIES

Due to its transparency, you can visualize the bath's inside. It improves temperature stability, and restricts evaporation and pollution

## MODELS

| Part No. | Temperature maximum °C | Capacity litres | Height / Width / Depth (usable) cm | Height / Width / Depth (exterior) cm | Power W | Weight Kg |
|----------|------------------------|-----------------|------------------------------------|--------------------------------------|---------|-----------|
| 6001120  | 99,9                   | 5               | 15 30 14                           | 20 37 27                             | 600     | 6         |
| 6001220  | 200                    | 5               | 15 30 14                           | 20 37 27                             | 1100    | 6         |
| 6001121  | 99,9                   | 12              | 15 30 30                           | 20 40 42                             | 1100    | 8         |
| 6001221  | 200                    | 12              | 15 30 30                           | 20 40 42                             | 1400    | 8         |
| 6001122  | 99,9                   | 20              | 15 48 30                           | 20 58 42                             | 1500    | 10        |
| 6001222  | 200                    | 20              | 15 48 30                           | 20 58 42                             | 2100    | 10        |
| 6001123  | 99,9                   | 45              | 15 62 50                           | 20 69 62                             | 2400    | 18        |
| 6001223  | 200                    | 45              | 15 62 50                           | 20 69 62                             | 3800    | 18        |

All baths with heating shelf and drainage. Faucet included.

## MODELS WITH MAKROLON® FOLDING LID

| Part No. | Temperature maximum °C | Capacity litres | Height / Width / Depth (usable) cm | Height / Width / Depth (exterior+lid) cm | Power W | Weight Kg |
|----------|------------------------|-----------------|------------------------------------|------------------------------------------|---------|-----------|
| 6001131  | 99,9                   | 12              | 15 30 30                           | 40 40 42                                 | 1100    | 8         |
| 6001132  | 99,9                   | 20              | 15 48 30                           | 40 58 42                                 | 1500    | 10        |
| 6001133  | 99,9                   | 45              | 15 62 50                           | 40 69 62                                 | 2400    | 18        |

All baths with heating shelf and drainage. Faucet included.

**ACCESSORIES.** See page 89.